

Content Determination of Potassium Dihydrogen Phosphate

Method for the alkalimetric content determination of Potassium Dihydrogen Phosphate

Potassium dihydrogen phosphate is a saline laxative, which is used to treat hyperpotassaemia from poor or malnourishment, or cases of high loss of potassium from vomiting or diarrhoea.

Sample	Potassium dihydrogen phosphate	Structure	
Substance	Potassium dihydrogen phosphate M = 136.1 g/mol	Literature	Ph. Eur. 1997; Monograph "Potassium dihydrogen phosphate"
Chemicals	Deionised water R	Sample preparation Weigh a sample of 0.24 to 0.28 g potassium dihydrogen phosphate, accurately to 0.1 mg, into a 100 mL titration beaker and dissolve in 40 mL deionised water.	
Titrant	c(NaOH) = 0.5 mol/l		
Standard	Potassium hydrogen phthalate (M521)		
Instruments	METTLER TOLEDO DL70ES/DL77 METTLER TOLEDO DL58/55/53/50Graphix		
Accessories	Titration beaker ME-101974 SP250 Peristaltic pump Rondo60 Sample changer		
Indication	DG111-SC	Remarks The previously determined loss on drying content is saved as the auxiliary value H2O.	
Chemistry	Deprotonation of the dihydrogen phosphate anions to monohydrogen phosphate $\text{H}_2\text{PO}_4^- + \text{OH}^- \rightarrow \text{HPO}_4^{2-} + \text{H}_2\text{O}$		
Calculation	$\text{R}(\%) = (\text{Q} \cdot \text{M} / (10 \cdot \text{m})) + \text{H2O}$ Q = Titrant consumption in mmol M = Molar mass K-dihydrogen phosphate m = Sample size in g H2O = % Loss on drying of KH_2PO_4	Waste disposal	Halogen-free organic solvents waste.
		Author	L. Pludra

Results

METTLER DL77

Titration V3.1

Dr. Kade Pharm. Fabrik GmbH

Tante Elfriede

Quality Control Laboratory

Method: G056 Cont. KH₂PO₄ Date: 14-12-1999 13:40
User: **Pludra** Measured: 14-12-1999 13:43

RESULTS

No.	ID1	Sample size and results
1/1	992222	0.2564 g R1 = 100.403 % KH ₂ PO ₄
1/2	992222	0.2453 g R1 = 100.270 % KH ₂ PO ₄
1/3	992222	0.2443 g R1 = 100.278 % KH ₂ PO ₄
1/4	992222	0.2581 g R1 = 100.294 % KH ₂ PO ₄
1/5	992222	0.2471 g R1 = 100.243 % KH ₂ PO ₄
1/6	992222	0.2553 g R1 = 100.171 % KH ₂ PO ₄


STATISTICS

Number of results R1 n = 6
Mean value $\bar{x}$ = 100.276 % KH₂PO₄
Standard deviation s = 0.075597 % KH₂PO₄
Rel. standard deviation srel = 0.075 %

Table of measured values

EQP	ET	Volume mL	Signal mV	1st Deriv. mV/mL	Time min:s
EQP1	ET1	0.0000	108.1		0:02
		9.9820	-25.7	-13.4	0:27
		14.9730	-54.5	-5.8	1:10
	ET2	17.4700	-78.7	-9.7	1:21
		17.9700	-86.8	-16.1	1:28
		18.3250	-94.3	-21.1	1:35
		18.6310	-102.5	-27.0	1:42
		18.8740	-111.4	-36.3	1:49
		19.0490	-120.2	-50.4	1:56
		19.1730	-128.5	-67.2	2:03
		19.2680	-137.0	-89.2	2:11
		19.3390	-145.5	-119.3	2:20
		19.3920	-153.2	-145.3	2:28
		19.4380	-160.2	-152.2	2:36
		19.4880	-168.8	-173.6	2:43
		19.5290	-175.8	-170.7	2:51
		19.5760	-183.2	-156.4	2:58
		19.6320	-190.9	-137.5	3:06
		19.6990	-198.0	-105.5	3:13
		19.8070	-206.9	-83.0	3:20
		19.9390	-215.2	-63.1	3:27
		20.0000	-218.4	-51.6	3:33

Titration curve



Method

Method: G056 Cont. KH2PO4
Version: 14-12-1999 13:40

Title
Method ID G056
Title Cont. KH2PO4
Date/time 14-Dec-1999 13:29

Sample
Number samples 6
Titration stand ST20 1
Entry type Weight m
Lower limit [g] 0.24
Upper limit [g] 0.26
ID1
Molar mass M 136.1
Equivalent number z 1
Temperature sensor Manual

Stir
Speed [%] 50
Time [s] 10

Titration
Titrant NaOH
Concentration [mol/L] 0.1
Sensor DG111-SC
Unit of meas. mV
Titration mode EQP
Predispensing 1 % nom. content
Metered amount [%] 90.0
Nominal content 100.0
Conversion const. M/(10*z)
Titrant addition DYN
dE(set) [mV] 8.0
Limits dV Absolute
dV(min) [mL] 0.02
dV(max) [mL] 0.5
Measure mode EQU
dE [mV] 0.5
dt [s] 1.0
t(min) [s] 2.0
t(max) [s] 20.0
Threshold 50.0
Maximum volume [mL] 20.0
Evaluation procedure Standard
Steepest jump only Yes
Stop for reevaluation Yes
Condition neq=0

Calculation
Result name KH2PO4
Formula $R = (Q \cdot C / m) + H2O$
Constant $C = M / (10 \cdot z)$
Result unit %
Decimal places 3

Rinse
Auxiliary reagent Water
Volume [mL] 30.0

Record
Output unit Computer
Raw results last sample Yes
Table of values Yes
E - V curve Yes

Statistics
R1 (i=index) R1
Standard deviation s Yes
Rel. standard deviation srel Yes

Record
Output unit Printer
Short-form method Yes
All results Yes