

Content Determination of Sodium Monohydrogen Phosphate

Method for the alkalimetric content determination of Sodium Monohydrogen Phosphate

Sodium monohydrogen phosphate is a saline laxative which is used for emptying the bowels, alkalising urine, hypercalcaemia and in cases of hypophosphataemia.

Sample	Na-monohydrogen phosphate	Structure	
Substance	Na-monohydrogen phosphate M = 142.0 g/mol	Literature	Ph. Eur. 1997; Monograph "Sodium Monohydrogen Phosphate"
Chemicals	Deionised water R	Sample preparation Weigh a sample of 0.24 to 0.28 g sodium monohydrogen phosphate, accurately to 0.1 mg, into a 100 mL titration beaker and dissolve in 40 mL deionised water.	
Titrant	c(NaOH) = 0.5 mol/L c(HCl) = 0.5 mol/L		
Standard	Potassium hydrogen phthalate (M521), THAM (M522)		
Instruments	METTLER TOLEDO DL70ES/DL77		
Accessories	Titration beaker ME-101974 SP250 Peristaltic pump Rondo60 Sample changer		
Indication	DG111-SC	Remarks The previously determined loss on drying content is saved as the auxiliary value H2O. The curve of the back titration shows two equivalence points so that the content of sodium monohydrogen phosphate as well as sodium dihydrogen phosphate can be determined on a parallel basis.	
Chemistry	1. Protonation of the monohydrogen phosphate to phosphoric acid through hydrochloric acid $\text{HPO}_4^{2-} + 2\text{H}_3\text{O}^+ \rightarrow \text{H}_3\text{PO}_4 + 2\text{H}_2\text{O}$ 2. Neutralisation of the excess hydrochloric acid as well as deprotonation to secondary phosphate in the first equivalence point $\text{H}_3\text{PO}_4 + \text{OH}^- \rightarrow \text{H}_2\text{PO}_4^- + \text{H}_2\text{O}$ 3. Deprotonation to primary phosphate in the second equivalence point $\text{H}_2\text{PO}_4^- + \text{OH}^- \rightarrow \text{HPO}_4^{2-}$		
Calculation	$\text{R}(\%) = (\text{QDISP} - \text{Q}) \cdot \text{M} \cdot \text{F} / (10 \cdot \text{m})$ QDISP= Total amount of dispensed HCl in mmol Q = Titrant consumption in mmol F = 100/(100-H2O) Formula for the loss on drying of Na ₂ HPO ₄ M = Molar mass Na ₂ HPO ₄ m = Sample size in g	Waste disposal Halogen-free organic solvents waste.	
		Author L. Pludra	

Results

METTLER DL77

Titration V3.1

Dr. Kade Pharm. Fabrik GmbH

Tante Elfriede

Quality Control Laboratory

Method: G031 Cont. Na2HPO4 Date: 14-12-1999 13:20
User: **Pludra** Measured: 14-12-1999 15:11

RESULTS

No.	ID1	Sample size and results			
1/1	992223	0.9178 g	R1 = 100.38 %	Na2HPO4	
1/2	992223	0.9477 g	R1 = 99.96 %	Na2HPO4	
1/3	992223	1.0699 g	R1 = 99.72 %	Na2HPO4	
1/4	992223	1.0368 g	R1 = 99.87 %	Na2HPO4	
1/5	992223	1.0046 g	R1 = 100.05 %	Na2HPO4	
1/6	992223	0.9735 g	R1 = 99.76 %	Na2HPO4	
			R3 = 19.00 %	Drying loss	


STATISTICS

Number of results R1 n = 6
Mean value $\bar{x}$ = 99.96 % Na2HPO4
Standard deviation s = 0.238792 % Na2HPO4
Rel. standard deviation srel = 0.239 %

Table of measured values

EQP	ET	Volume mL	Signal mV	1st Deriv. mV/mL	Time min:s
EQP1	ET1	0.0000	320.3		0:03
		5.7140	301.0	-3.4	0:13
	ET2	8.5700	285.5	-5.4	0:20
		10.0000	275.7	-6.9	0:26
		10.2000	273.8	-9.5	0:29
		10.4000	272.1	-8.8	0:33
		10.6000	270.4	-8.4	0:36
		10.8000	268.5	-9.5	0:40
		11.0000	267.0	-7.4	0:44
		\\	\\	\\	\\
		14.5440	173.9	-148.1	2:10
		14.5880	165.2	-198.9	2:19
		14.6200	158.0	-225.3	2:25
		14.6500	150.7	-242.7	2:31
		14.6800	142.3	-280.0	2:38
		14.7040	136.8	-227.5	2:45
		14.7480	128.5	-189.3	2:51
		14.8000	119.5	-173.7	2:59
		\\	\\	\\	\\
		24.4040	-62.4	-44.8	6:53
	24.5400	-68.9	-48.4	6:57	
	24.6920	-80.1	-73.2	7:02	
	24.7720	-87.1	-87.5	7:08	
	EQP2		24.8500	-96.9	-126.5
		24.8980	-105.1	-170.6	7:20
		24.9340	-112.1	-192.5	7:26
		24.9700	-119.8	-213.9	7:33
		25.0040	-127.9	-238.8	7:39
		25.0340	-134.9	-233.3	7:45
		25.0680	-142.2	-216.2	7:52

Titration curve



Method

Method: G031 Cont. Na2HPO4
Version: 14-12-1999 13:20

Title
Method ID G031
Title Cont. Na2HPO4
Date/time 14-Dec-1999 13:20

Sample
Number samples 6
Titration stand ST20 1
Entry type Weight m
Lower limit [g] 0.9
Upper limit [g] 1.1
ID1
Molar mass M 100.0
Equivalent number z 1
Temperature sensor Manual

Dispense
Titrant HCl
Concentration [mol/L] 0.5
Volume [mL] 25

Stir
Speed [%] 50
Time [s] 10

Titration
Titrant NaOH
Concentration [mol/L] 0.5
Sensor DG111-SC
Unit of meas. mV
Titration mode EQP
Predispensing 1 mL
Volume [mL] 10
Titrant addition DYN
dE(set) [mV] 8.0
Limits dV Absolute
dV(min) [mL] 0.02
dV(max) [mL] 0.2
Measure mode EQU
dE [mV] 0.5
dt [s] 1.0
t(min) [s] 3.0
t(max) [s] 30.0
Threshold 100.0
Maximum volume [mL] 30.0
Termination after n EQPs Yes
n = 2
Evaluation procedure Standard
Stop for reevaluation Yes
Condition neq<2

Calculation
Result name Na2HPO4
Formula $R1=1420 \cdot (QDISP-Q1) / C$
Constant $C=m \cdot (100-H20)$
Result unit %
Decimal places 2

Calculation
Result name NaH2PO4
Formula $R2=(Q1+Q2-QDISP) / C2$
Constant $C2=(QDISP-Q1) / 100$
Result unit %
Decimal places 2
Condition Yes
Condition $Q1+Q2>QDISP$

Rinse
Auxiliary reagent Water
Volume [mL] 30.0

Record
Output unit Computer
Raw results last sample Yes
Table of values Yes
E - V curve Yes

Statistics
Ri (i=index) R1
Standard deviation s Yes
Rel. standard deviation srel Yes

Statistics
Ri (i=index) R2
Standard deviation s Yes
Rel. standard deviation srel Yes

Calculation
Result name Loss on drying
Formula $R3=H20$
Constant
Result unit %
Decimal places 2

Record
Output unit Printer
All results Yes