

**Method:** TitreCericSulphate    **Titer** 0.1mol/L CeSO4    **8/10/2012 10:44:00 AM**  
**Start time:** 8/10/2012 11:11:10 AM

### Sample data

No.	Comment / ID	Start time	Sample size	Corr. f	Density
1/6	Disodium oxalat	8/10/2012 11:11:11 AM	0.03293 g	1.0	0 g/mL
2/6	Disodium oxalat	8/10/2012 11:17:04 AM	0.03222 g	1.0	0 g/mL
3/6	Disodium oxalat	8/10/2012 11:22:43 AM	0.03566 g	1.0	0 g/mL
4/6	Disodium oxalat	8/10/2012 11:29:20 AM	0.03385 g	1.0	0 g/mL
5/6	Disodium oxalat	8/10/2012 11:35:29 AM	0.03482 g	1.0	0 g/mL
6/6	Disodium oxalat	8/10/2012 11:41:53 AM	0.03653 g	1.0	0 g/mL

### Results

No.	Comment / ID	Start time	Sample size and results			
1/6	Disodium oxalat	8/10/2012 11:11:11 AM	0.03293	g		
			R1 = 0.99782	--		Titer
2/6	Disodium oxalat	8/10/2012 11:17:04 AM	0.03222	g		
			R1 = 0.99949	--		Titer
3/6	Disodium oxalat	8/10/2012 11:22:43 AM	0.03566	g		
			R1 = 0.99963	--		Titer
4/6	Disodium oxalat	8/10/2012 11:29:20 AM	0.03385	g		
			R1 = 0.99755	--		Titer
5/6	Disodium oxalat	8/10/2012 11:35:29 AM	0.03482	g		
			R1 = 0.99690	--		Titer
6/6	Disodium oxalat	8/10/2012 11:41:53 AM	0.03653	g		
			R1 = 0.99707	--		Titer
-/-			R2 = 0.9981	--		Mean Titer

### Titer

Titer      0.99808

### Series comment

### Statistics

Rx	Name	n	Mean value	Unit	s	srel [%]
R1	Titer	6	0.99808	--	0.00120	0.120
R2	Mean Titer	1	0.9981	--	NaN	NaN

### Raw data

### Sample

No.      1/6  
Standard      Disodium oxalat  
Type of standard      solid  
Comment  
Titration stand      Rondo60/1A  
Weight      m = 0.03293 g  
Correction factor      f = 1.0

**Method:** TitreCericSulphate Titer 0.1mol/L CeSO4  
**Start time:** 8/10/2012 11:11:10 AM 8/10/2012 10:44:00 AM

Purity p = 100.00 %  
Temperature T = 25.0 oC  
Sample start 8/10/2012 11:11:11 AM  
Sample end 8/10/2012 11:17:04 AM

**Dispense (normal) [1]**

Titrant H2SO4 cDi = 1 mol/L TITERDi = 1.0  
Disp. volume VENDDi = 5.0 mL  
Disp. amount QENDDi = 5.00 mmol  
Time 0:13 min

**Measure (normal) [1]**

Sensor DH 100  
Temperature DH 100 241.5 oC  
Measured value  
Time tMe = 0:01 min

**EQP titration [1]**

Titrant Ce(SO4)2 c = 0.1 mol/L TITER = 0.99088  
Sensor DM140-SC  
Start potential EST = 603.2 mV  
Predispense EPD = 596.2 mV  
VPD = 4.0000 mL  
nEQ = 1  
No. of EQPs and cand. EQP1  
Consumption VEQ1 = 4.925293 mL  
Q1 = 0.488037 mmol  
EEQ1 = 755.4 mV  
EHN1 = 625.2 mV  
Excess VEX = 0.474707 mL  
QEX = 0.047038 mmol  
End VEND = 5.4000 mL  
QEND = 0.535075 mmol  
Termination at EQPs  
Time t = 3:17 min

**Calculation**

Result R1 = 0.99782 -- Titer  
Formula  $R1 = m / (VEQ * c * C)$   
Constant  $M / (10 * p * z)$   
C = 0.067005  
Molar mass M[Disodium oxalat] = 134.01 g/mol  
Equivalent number z[Disodium oxalat] = 2  
Duration tUSE = 05:19 min

**Measured values EQP titration [1]**

Titrant Ce(SO4)2 c = 0.1 mol/L TITER = 0.99088  
Sensor DM140-SC  
Sample 1/6


Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature oC
0.0000	NaN	603.2	NaN	NaN	0	25.0
4.0000	4.0000	639.0	35.8	NaN	12	25.0
4.0500	0.0500	596.6	-42.4	NaN	48	25.0
4.1000	0.0500	596.1	-0.5	NaN	52	25.0
4.1500	0.0500	595.4	-0.7	NaN	58	25.0
4.2000	0.0500	594.8	-0.6	-73.30	63	25.0
4.2500	0.0500	594.6	-0.2	44.63	68	25.0
4.3000	0.0500	594.4	-0.2	1.30	73	25.0
4.3500	0.0500	594.5	0.1	6.40	78	25.0
4.4000	0.0500	595.3	0.8	10.54	83	25.0
4.4500	0.0500	595.9	0.6	13.63	88	25.0

**Method:** TitreCericSulphate **Titer** 0.1mol/L CeSO4  
**Start time:** 8/10/2012 11:11:10 AM

8/10/2012 10:44:00 AM

	Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature oC
	4.5000	0.0500	596.7	0.8	20.72	93	25.0
	4.5500	0.0500	597.7	1.0	26.12	98	25.0
	4.6000	0.0500	599.3	1.6	28.13	103	25.0
	4.6500	0.0500	601.5	2.2	6.71	108	25.0
	4.7000	0.0500	605.3	3.8	-69.84	113	25.0
	4.7500	0.0500	608.4	3.1	-14.44	118	25.0
	4.8000	0.0500	615.3	6.9	664.92	123	25.0
	4.8500	0.0500	630.3	15.0	1529.96	128	25.0
	4.9000	0.0500	680.1	49.8	1946.35	133	25.0
EQP1	4.925293	NaN	755.4	NaN	1951.10	NaN	NaN
	4.9500	0.0500	829.0	148.9	1723.07	138	25.0
	5.0000	0.0500	1037.2	208.2	976.40	143	25.0
	5.0500	0.0500	1006.8	-30.4	-4.61	148	25.0
	5.1000	0.0500	883.5	-123.3	-868.00	153	25.0
	5.1500	0.0500	814.1	-69.4	-1254.08	158	25.0
	5.2000	0.0500	787.7	-26.4	NaN	163	25.0
	5.2500	0.0500	782.0	-5.7	NaN	168	25.0
	5.3000	0.0500	782.1	0.1	NaN	173	25.0
	5.3500	0.0500	783.2	1.1	NaN	178	25.0
	5.4000	0.0500	787.4	4.2	NaN	183	25.0


**E - V curve** EQP titration [1]  
Sample 1/6



**Method:** TitreCericSulphate Titer 0.1mol/L CeSO<sub>4</sub>  
**Start time:** 8/10/2012 11:11:10 AM

8/10/2012 10:44:00 AM

**dE/dV - V curve** EQP titration [1]  
Sample 1/6



## Raw data

### Sample

No.	2/6
Standard	Disodium oxalat
Type of standard	solid
Comment	
Titration stand	Rondo60/1A
Weight	m = 0.03222 g
Correction factor	f = 1.0
Purity	p = 100.00 %
Temperature	T = 25.0 oC
Sample start	8/10/2012 11:17:04 AM
Sample end	8/10/2012 11:22:43 AM

### Dispense (normal) [1]

Titrant	H <sub>2</sub> SO <sub>4</sub> cDi = 1 mol/L TITERDi = 1.0
Disp. volume	VENDDi = 5.0 mL
Disp. amount	QENDDi = 5.00 mmol
Time	0:12 min

### Measure (normal) [1]

Sensor	DH 100
Temperature	DH 100 241.5 oC
Measured value	
Time	tMe = 0:01 min

### EQP titration [1]

Titrant	Ce(SO <sub>4</sub> ) <sub>2</sub> c = 0.1 mol/L TITER = 0.99088
Sensor	DM140-SC
Start potential	EST = 598.1 mV
Predispense	EPD = 597.4 mV
	VPD = 4.0000 mL
	nEQ = 1
No. of EQPs and cand.	
Consumption	VEQ1 = 4.811039 mL

**Method:** TitreCericSulphate **Titer** 0.1mol/L CeSO4 **8/10/2012 10:44:00 AM**  
**Start time:** 8/10/2012 11:11:10 AM

Excess  
End  
Termination at Time

EQPs  
t = 3:01 min

Q1 = 0.476716 mmol  
EEQ1 = 738.4 mV  
EHN1 = 621.6 mV  
VEX = 0.438961 mL  
QEX = 0.043496 mmol  
VEND = 5.2500 mL  
QEND = 0.520212 mmol

# **Calculation**

Result  
Formula  
Constant  
Molar mass  
Equivalent number  
Duration

R1 = 0.99949 -- Titer  
R1=m/(VEQ\*c\*C)  
M/(10\*p\*z)  
C = 0.067005  
M[Disodium oxalat] = 134.01 g/mol  
z[Disodium oxalat] = 2  
tUSE = 05:04 min

# **Measured values EQP titration [1]**

Titrant  
Sensor  
Sample


Ce(SO4)2 c = 0.1 mol/L TITER = 0.99088  
DM140-SC  
2/6

	Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature oC
	0.0000	NaN	598.1	NaN	NaN	0	25.0
	4.0000	4.0000	637.1	39.0	NaN	12	25.0
	4.0500	0.0500	597.2	-39.9	NaN	48	25.0
	4.1000	0.0500	596.3	-0.9	NaN	53	25.0
	4.1500	0.0500	595.8	-0.5	NaN	58	25.0
	4.2000	0.0500	595.6	-0.2	-65.66	63	25.0
	4.2500	0.0500	595.7	0.1	47.83	68	25.0
	4.3000	0.0500	595.8	0.1	7.08	73	25.0
	4.3500	0.0500	596.5	0.7	11.27	78	25.0
	4.4000	0.0500	596.8	0.3	15.32	83	25.0
	4.4500	0.0500	598.3	1.5	19.48	88	25.0
	4.5000	0.0500	599.4	1.1	18.29	93	25.0
	4.5500	0.0500	601.4	2.0	1.88	98	25.0
	4.6000	0.0500	604.2	2.8	-101.26	103	25.0
	4.6500	0.0500	609.0	4.8	165.00	108	25.0
	4.7000	0.0500	617.6	8.6	924.34	113	25.0
	4.7500	0.0500	639.0	21.4	1500.71	118	25.0
	4.8000	0.0500	695.4	56.4	1597.46	123	25.0
EQP1	4.811039	NaN	738.4	NaN	1618.99	NaN	NaN
	4.8500	0.0500	890.2	194.8	1209.44	128	25.0
	4.9000	0.0500	987.8	97.6	510.61	133	25.0
	4.9500	0.0500	880.5	-107.3	-256.32	138	25.0
	5.0000	0.0500	813.8	-66.7	-806.26	143	25.0
	5.0500	0.0500	786.2	-27.6	NaN	148	25.0
	5.1000	0.0500	781.5	-4.7	NaN	153	25.0
	5.1500	0.0500	781.5	0.0	NaN	158	25.0
	5.2000	0.0500	783.8	2.3	NaN	163	25.0
	5.2500	0.0500	785.7	1.9	NaN	168	25.0


**Method:** TitreCericSulphate Titer 0.1mol/L CeSO<sub>4</sub>  
**Start time:** 8/10/2012 11:11:10 AM

8/10/2012 10:44:00 AM

**E - V curve** EQP titration [1]  
Sample 2/6



**dE/dV - V curve** EQP titration [1]  
Sample 2/6



## Raw data

### Sample

No.	3/6
Standard	Disodium oxalat
Type of standard	solid
Comment	
Titration stand	Rondo60/1A
Weight	m = 0.03566 g
Correction factor	f = 1.0
Purity	p = 100.00 %

**Method:** TitreCericSulphate Titer 0.1mol/L CeSO4  
**Start time:** 8/10/2012 11:11:10 AM 8/10/2012 10:44:00 AM

Temperature T = 25.0 oC  
Sample start 8/10/2012 11:22:43 AM  
Sample end 8/10/2012 11:29:19 AM

**Dispense (normal) [1]**

Titrant H2SO4 cDi = 1 mol/L TITERDi = 1.0  
Disp. volume VENDDi = 5.0 mL  
Disp. amount QENDDi = 5.00 mmol  
Time 0:12 min

**Measure (normal) [1]**

Sensor DH 100  
Temperature DH 100 241.5 oC  
Measured value  
Time tMe = 0:01 min

**EQP titration [1]**

Titrant Ce(SO4)2 c = 0.1 mol/L TITER = 0.99088  
Sensor DM140-SC  
Start potential EST = 595.2 mV  
Predispense EPD = 590.3 mV  
VPD = 4.0000 mL  
nEQ = 1  
No. of EQPs and cand. EQP1 VEQ1 = 5.323951 mL  
Consumption Q1 = 0.527540 mmol  
EEQ1 = 732.9 mV  
EHN1 = 616.0 mV  
VEX = 0.476049 mL  
QEX = 0.047171 mmol  
VEND = 5.8000 mL  
QEND = 0.574710 mmol  
Excess  
End  
Termination at EQPs  
Time t = 3:58 min

**Calculation**

Result R1 = 0.99963 -- Titer  
Formula  $R1 = m / (VEQ * c * C)$   
Constant  $M / (10 * p * z)$   
C = 0.067005  
Molar mass M[Disodium oxalat] = 134.01 g/mol  
Equivalent number z[Disodium oxalat] = 2  
Duration tUSE = 06:01 min

**Measured values EQP titration [1]**

Titrant Ce(SO4)2 c = 0.1 mol/L TITER = 0.99088  
Sensor DM140-SC  
Sample 3/6


Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature oC
0.0000	NaN	595.2	NaN	NaN	0	25.0
4.0000	4.0000	626.5	31.3	NaN	12	25.0
4.0500	0.0500	589.7	-36.8	NaN	48	25.0
4.1000	0.0500	588.1	-1.6	NaN	53	25.0
4.1500	0.0500	587.2	-0.9	NaN	58	25.0
4.2000	0.0500	586.1	-1.1	-71.54	63	25.0
4.2500	0.0500	585.8	-0.3	33.49	68	25.0
4.3000	0.0500	585.5	-0.3	-3.92	73	25.0
4.3500	0.0500	585.5	0.0	-1.02	78	25.0
4.4000	0.0500	585.5	0.0	1.91	83	25.0
4.4500	0.0500	585.5	0.0	3.31	88	25.0
4.5000	0.0500	585.6	0.1	5.55	93	25.0

**Method:** TitreCericSulphate **Titer** 0.1mol/L CeSO4  
**Start time:** 8/10/2012 11:11:10 AM

8/10/2012 10:44:00 AM

	Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature oC
	4.5500	0.0500	586.3	0.7	8.79	98	25.0
	4.6000	0.0500	586.7	0.4	12.87	103	25.0
	4.6500	0.0500	587.2	0.5	16.42	108	25.0
	4.7000	0.0500	588.2	1.0	19.36	113	25.0
	4.7500	0.0500	589.5	1.3	22.39	118	25.0
	4.8000	0.0500	590.6	1.1	25.00	123	25.0
	4.8500	0.0500	591.8	1.2	25.89	128	25.0
	4.9000	0.0500	593.3	1.5	28.86	133	25.0
	4.9500	0.0500	595.0	1.7	33.78	138	25.0
	5.0000	0.0500	596.7	1.7	42.23	143	25.0
	5.0500	0.0500	599.3	2.6	24.00	148	25.0
	5.1000	0.0500	603.3	4.0	-33.57	153	25.0
	5.1500	0.0500	608.0	4.7	3.28	158	25.0
	5.2000	0.0500	615.8	7.8	619.27	163	25.0
	5.2500	0.0500	627.7	11.9	1332.38	168	25.0
	5.3000	0.0500	673.4	45.7	1628.50	173	25.0
EQP1	5.323951	NaN	732.9	NaN	1633.83	NaN	NaN
	5.3500	0.0500	797.7	124.3	1409.34	178	25.0
	5.4000	0.0500	985.1	187.4	792.27	183	25.0
	5.4500	0.0500	933.1	-52.0	5.83	188	25.0
	5.5000	0.0500	832.6	-100.5	-681.83	193	25.0
	5.5500	0.0500	790.2	-42.4	-955.85	198	25.0
	5.6000	0.0500	776.7	-13.5	NaN	204	25.0
	5.6500	0.0500	774.6	-2.1	NaN	208	25.0
	5.7000	0.0500	774.9	0.3	NaN	214	25.0
	5.7500	0.0500	780.8	5.9	NaN	218	25.0
	5.8000	0.0500	783.3	2.5	NaN	224	25.0

**E - V curve** EQP titration [1]  
Sample 3/6






**Method:** TitreCericSulphate Titer 0.1mol/L CeSO<sub>4</sub>  
**Start time:** 8/10/2012 11:11:10 AM

8/10/2012 10:44:00 AM

**dE/dV - V curve EQP titration [1]**  
Sample 3/6



## Raw data

### Sample

No.	4/6
Standard	Disodium oxalat
Type of standard	solid
Comment	
Titration stand	Rondo60/1A
Weight	m = 0.03385 g
Correction factor	f = 1.0
Purity	p = 100.00 %
Temperature	T = 25.0 oC
Sample start	8/10/2012 11:29:20 AM
Sample end	8/10/2012 11:35:29 AM

### Dispense (normal) [1]

Titrant	H <sub>2</sub> SO <sub>4</sub> cDi = 1 mol/L TITERDi = 1.0
Disp. volume	VENDDi = 5.0 mL
Disp. amount	QENDDi = 5.00 mmol
Time	0:12 min

### Measure (normal) [1]

Sensor	DH 100
Temperature	DH 100 241.5 oC
Measured value	
Time	tMe = 0:01 min

### EQP titration [1]

Titrant	Ce(SO <sub>4</sub> ) <sub>2</sub> c = 0.1 mol/L TITER = 0.99088
Sensor	DM140-SC
Start potential	EST = 598.6 mV
Predispense	EPD = 595.7 mV
	VPD = 4.0000 mL
	nEQ = 1
No. of EQPs and cand.	
Consumption	VEQ1 = 5.064252 mL

**Method:** TitreCericSulphate **Titer** 0.1mol/L CeSO4 **8/10/2012 10:44:00 AM**  
**Start time:** 8/10/2012 11:11:10 AM

Excess  
End  
Termination at Time

EQPs  
t = 3:28 min

Q1 = 0.501807 mmol  
EEQ1 = 734.2 mV  
EHN1 = 618.2 mV  
VEX = 0.435748 mL  
QEX = 0.043177 mmol  
VEND = 5.5000 mL  
QEND = 0.544984 mmol

# **Calculation**

Result  
Formula  
Constant  
Molar mass  
Equivalent number  
Duration

R1 = 0.99755 -- Titer  
R1=m/(VEQ\*c\*C)  
M/(10\*p\*z)  
C = 0.067005  
M[Disodium oxalat] = 134.01 g/mol  
z[Disodium oxalat] = 2  
tUSE = 05:32 min

# **Measured values EQP titration [1]**

Titrant  
Sensor  
Sample


Ce(SO4)2 c = 0.1 mol/L TITER = 0.99088  
DM140-SC  
4/6

	Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature oC
	0.0000	NaN	598.6	NaN	NaN	0	25.0
	4.0000	4.0000	629.6	31.0	NaN	12	25.0
	4.0500	0.0500	595.0	-34.6	NaN	49	25.0
	4.1000	0.0500	594.0	-1.0	NaN	54	25.0
	4.1500	0.0500	593.2	-0.8	NaN	59	25.0
	4.2000	0.0500	592.6	-0.6	-64.06	64	25.0
	4.2500	0.0500	592.2	-0.4	32.84	69	25.0
	4.3000	0.0500	592.0	-0.2	-2.22	74	25.0
	4.3500	0.0500	592.0	0.0	2.54	79	25.0
	4.4000	0.0500	592.1	0.1	6.88	84	25.0
	4.4500	0.0500	592.7	0.6	9.87	89	25.0
	4.5000	0.0500	593.2	0.5	12.66	94	25.0
	4.5500	0.0500	594.2	1.0	15.60	99	25.0
	4.6000	0.0500	594.8	0.6	20.30	104	25.0
	4.6500	0.0500	595.7	0.9	22.62	109	25.0
	4.7000	0.0500	597.3	1.6	27.89	114	25.0
	4.7500	0.0500	599.1	1.8	27.24	119	25.0
	4.8000	0.0500	601.4	2.3	-2.57	124	25.0
	4.8500	0.0500	603.4	2.0	-86.89	129	25.0
	4.9000	0.0500	608.4	5.0	76.62	134	25.0
	4.9500	0.0500	614.4	6.0	818.91	139	25.0
	5.0000	0.0500	630.7	16.3	1451.21	144	25.0
	5.0500	0.0500	687.6	56.9	1604.78	149	25.0
EQP1	5.064252	NaN	734.2	NaN	1619.77	NaN	NaN
	5.1000	0.0500	851.2	163.6	1263.73	154	25.0
	5.1500	0.0500	995.7	144.5	589.94	159	25.0
	5.2000	0.0500	889.5	-106.2	-180.44	164	25.0
	5.2500	0.0500	811.3	-78.2	-767.82	169	25.0
	5.3000	0.0500	783.6	-27.7	NaN	174	25.0
	5.3500	0.0500	777.9	-5.7	NaN	179	25.0
	5.4000	0.0500	778.8	0.9	NaN	184	25.0
	5.4500	0.0500	779.2	0.4	NaN	189	25.0
	5.5000	0.0500	781.2	2.0	NaN	194	25.0


**Method:** TitreCericSulphate Titer 0.1mol/L CeSO<sub>4</sub>  
**Start time:** 8/10/2012 11:11:10 AM

8/10/2012 10:44:00 AM

**E - V curve** EQP titration [1]  
Sample 4/6



**dE/dV - V curve** EQP titration [1]  
Sample 4/6



## Raw data

### Sample

No.	5/6
Standard	Disodium oxalat
Type of standard	solid
Comment	
Titration stand	Rondo60/1A
Weight	m = 0.03482 g
Correction factor	f = 1.0
Purity	p = 100.00 %

**Method:** TitreCericSulphate Titer 0.1mol/L CeSO4 8/10/2012 10:44:00 AM  
**Start time:** 8/10/2012 11:11:10 AM

Temperature T = 25.0 oC  
Sample start 8/10/2012 11:35:29 AM  
Sample end 8/10/2012 11:41:53 AM

**Dispense (normal) [1]**

Titrant H2SO4 cDi = 1 mol/L TITERDi = 1.0  
Disp. volume VENDDi = 5.0 mL  
Disp. amount QENDDi = 5.00 mmol  
Time 0:13 min

**Measure (normal) [1]**

Sensor DH 100  
Temperature DH 100 241.5 oC  
Measured value  
Time tMe = 0:01 min

**EQP titration [1]**

Titrant Ce(SO4)2 c = 0.1 mol/L TITER = 0.99088  
Sensor DM140-SC  
Start potential EST = 599.6 mV  
Predispense EPD = 594.8 mV  
VPD = 4.0000 mL  
nEQ = 1  
No. of EQPs and cand. EQP1 VEQ1 = 5.212798 mL  
Consumption Q1 = 0.516526 mmol  
EEQ1 = 738.7 mV  
EHN1 = 620.1 mV  
Excess VEX = 0.437202 mL  
QEX = 0.043321 mmol  
End VEND = 5.6500 mL  
QEND = 0.559847 mmol  
Termination at EQPs  
Time t = 3:43 min

**Calculation**

Result R1 = 0.99690 -- Titer  
Formula  $R1 = m / (VEQ * c * C)$   
Constant  $M / (10 * p * z)$   
C = 0.067005  
Molar mass M[Disodium oxalat] = 134.01 g/mol  
Equivalent number z[Disodium oxalat] = 2  
Duration tUSE = 05:45 min

**Measured values EQP titration [1]**

Titrant Ce(SO4)2 c = 0.1 mol/L TITER = 0.99088  
Sensor DM140-SC  
Sample 5/6


Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature oC
0.0000	NaN	599.6	NaN	NaN	0	25.0
4.0000	4.0000	631.0	31.4	NaN	12	25.0
4.0500	0.0500	594.3	-36.7	NaN	48	25.0
4.1000	0.0500	593.2	-1.1	NaN	53	25.0
4.1500	0.0500	592.2	-1.0	NaN	58	25.0
4.2000	0.0500	591.4	-0.8	-68.96	63	25.0
4.2500	0.0500	591.1	-0.3	34.84	68	25.0
4.3000	0.0500	590.9	-0.2	-2.66	73	25.0
4.3500	0.0500	590.9	0.0	0.35	78	25.0
4.4000	0.0500	591.0	0.1	3.27	83	25.0
4.4500	0.0500	591.1	0.1	5.55	88	25.0
4.5000	0.0500	591.4	0.3	8.70	93	25.0

**Method:** TitreCericSulphate **Titer** 0.1mol/L CeSO4  
**Start time:** 8/10/2012 11:11:10 AM

8/10/2012 10:44:00 AM

	Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature oC
	4.5500	0.0500	592.0	0.6	11.25	98	25.0
	4.6000	0.0500	592.8	0.8	15.68	103	25.0
	4.6500	0.0500	593.4	0.6	17.48	108	25.0
	4.7000	0.0500	594.5	1.1	20.20	113	25.0
	4.7500	0.0500	595.5	1.0	23.25	118	25.0
	4.8000	0.0500	597.2	1.7	26.56	123	25.0
	4.8500	0.0500	597.9	0.7	32.61	128	25.0
	4.9000	0.0500	600.5	2.6	34.06	133	25.0
	4.9500	0.0500	602.6	2.1	-6.08	138	25.0
	5.0000	0.0500	605.7	3.1	-87.83	143	25.0
	5.0500	0.0500	610.6	4.9	87.54	148	25.0
	5.1000	0.0500	616.7	6.1	864.99	153	25.0
	5.1500	0.0500	631.6	14.9	1501.56	158	25.0
	5.2000	0.0500	696.1	64.5	1641.49	163	25.0
EQP1	5.212798	NaN	738.7	NaN	1659.89	NaN	NaN
	5.2500	0.0500	862.6	166.5	1271.95	168	25.0
	5.3000	0.0500	1010.5	147.9	564.30	173	25.0
	5.3500	0.0500	890.9	-119.6	-236.01	178	25.0
	5.4000	0.0500	813.4	-77.5	-836.18	183	25.0
	5.4500	0.0500	783.0	-30.4	NaN	188	25.0
	5.5000	0.0500	777.1	-5.9	NaN	193	25.0
	5.5500	0.0500	776.7	-0.4	NaN	198	25.0
	5.6000	0.0500	777.1	0.4	NaN	203	25.0
	5.6500	0.0500	778.3	1.2	NaN	208	25.0


**E - V curve** **EQP titration [1]**  
Sample 5/6



**Method:** TitreCericSulphate Titer 0.1mol/L CeSO<sub>4</sub>  
**Start time:** 8/10/2012 11:11:10 AM

8/10/2012 10:44:00 AM

**dE/dV - V curve** EQP titration [1]  
Sample 5/6



## Raw data

### Sample

No. 6/6  
Standard Disodium oxalat  
Type of standard solid  
Comment  
Titration stand Rondo60/1A  
Weight m = 0.03653 g  
Correction factor f = 1.0  
Purity p = 100.00 %  
Temperature T = 25.0 oC  
Sample start 8/10/2012 11:41:53 AM  
Sample end 8/10/2012 11:48:44 AM

### Dispense (normal) [1]

Titrant H<sub>2</sub>SO<sub>4</sub> cDi = 1 mol/L TITERDi = 1.0  
Disp. volume VENDDi = 5.0 mL  
Disp. amount QENDDi = 5.00 mmol  
Time 0:12 min

### Measure (normal) [1]

Sensor DH 100  
Temperature DH 100 241.5 oC  
Measured value  
Time tMe = 0:01 min

### EQP titration [1]

Titrant Ce(SO<sub>4</sub>)<sub>2</sub> c = 0.1 mol/L TITER = 0.99088  
Sensor DM140-SC  
Start potential EST = 596.1 mV  
Predispense EPD = 592.2 mV  
VPD = 4.0000 mL  
No. of EQPs and cand. nEQ = 1  
Consumption EQP1 VEQ1 = 5.467850 mL

**Method:** TitreCericSulphate **Titer** 0.1mol/L CeSO4 **8/10/2012 10:44:00 AM**  
**Start time:** 8/10/2012 11:11:10 AM

Excess  
End  
Termination at Time

EQPs  
t = 4:09 min

Q1 = 0.541798 mmol  
EEQ1 = 727.8 mV  
EHN1 = 611.5 mV  
VEX = 0.432150 mL  
QEX = 0.042821 mmol  
VEND = 5.9000 mL  
QEND = 0.584619 mmol

# **Calculation**

Result  
Formula  
Constant  
Molar mass  
Equivalent number  
Duration

R1 = 0.99707 -- Titer  
R1=m/(VEQ\*c\*C)  
M/(10\*p\*z)  
C = 0.067005  
M[Disodium oxalat] = 134.01 g/mol  
z[Disodium oxalat] = 2  
tUSE = 06:12 min

# **Measured values EQP titration [1]**

Titrant  
Sensor  
Sample


Ce(SO4)2 c = 0.1 mol/L TITER = 0.99088  
DM140-SC  
6/6

	Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature oC
	0.0000	NaN	596.1	NaN	NaN	0	25.0
	4.0000	4.0000	618.7	22.6	NaN	12	25.0
	4.0500	0.0500	591.8	-26.9	NaN	48	25.0
	4.1000	0.0500	590.8	-1.0	NaN	53	25.0
	4.1500	0.0500	589.9	-0.9	NaN	58	25.0
	4.2000	0.0500	589.1	-0.8	-55.05	63	25.0
	4.2500	0.0500	588.5	-0.6	19.32	68	25.0
	4.3000	0.0500	588.0	-0.5	-7.53	73	25.0
	4.3500	0.0500	587.8	-0.2	-4.20	78	25.0
	4.4000	0.0500	587.7	-0.1	-0.48	83	25.0
	4.4500	0.0500	587.6	-0.1	3.05	88	25.0
	4.5000	0.0500	587.9	0.3	6.53	93	25.0
	4.5500	0.0500	588.3	0.4	8.78	98	25.0
	4.6000	0.0500	588.9	0.6	10.14	103	25.0
	4.6500	0.0500	589.4	0.5	11.06	108	25.0
	4.7000	0.0500	590.0	0.6	11.35	113	25.0
	4.7500	0.0500	590.4	0.4	12.76	118	25.0
	4.8000	0.0500	591.1	0.7	14.73	123	25.0
	4.8500	0.0500	592.0	0.9	15.77	128	25.0
	4.9000	0.0500	592.8	0.8	16.35	133	25.0
	4.9500	0.0500	593.9	1.1	20.17	138	25.0
	5.0000	0.0500	594.8	0.9	24.62	143	25.0
	5.0500	0.0500	595.7	0.9	29.29	148	25.0
	5.1000	0.0500	597.6	1.9	31.87	153	25.0
	5.1500	0.0500	600.3	2.7	34.70	158	25.0
	5.2000	0.0500	602.1	1.8	18.55	163	25.0
	5.2500	0.0500	604.6	2.5	-100.31	168	25.0
	5.3000	0.0500	608.3	3.7	38.90	173	25.0
	5.3500	0.0500	615.6	7.3	715.05	178	25.0
	5.4000	0.0500	627.6	12.0	1337.38	183	25.0
	5.4500	0.0500	667.6	40.0	1531.13	188	25.0
EQP1	5.467850	NaN	727.8	NaN	1540.46	NaN	NaN
	5.5000	0.0500	836.3	168.7	1252.63	193	25.0
	5.5500	0.0500	969.2	132.9	641.46	198	25.0
	5.6000	0.0500	885.6	-83.6	-87.06	203	25.0
	5.6500	0.0500	810.8	-74.8	-677.98	208	25.0
	5.7000	0.0500	781.0	-29.8	NaN	213	25.0
	5.7500	0.0500	775.4	-5.6	NaN	219	25.0
	5.8000	0.0500	775.3	-0.1	NaN	224	25.0
	5.8500	0.0500	777.2	1.9	NaN	228	25.0
	5.9000	0.0500	780.4	3.2	NaN	234	25.0


**Method:** TitreCericSulphate Titer 0.1mol/L CeSO<sub>4</sub>  
**Start time:** 8/10/2012 11:11:10 AM

8/10/2012 10:44:00 AM

**E - V curve** EQP titration [1]  
Sample 6/6



**dE/dV - V curve** EQP titration [1]  
Sample 6/6



## Raw data

### Calculation

Result	R <sup>2</sup> = 0.9981 -- Mean Titer
Formula	R <sup>2</sup> = Mean[R <sup>1</sup> ]
Constant	1
	C = 1
Molar mass	M[None] = 1 g/mol
Equivalent number	z[None] = 1



---

<b>Method:</b>	<b>TitreCericSulphate</b>	<b>Titer 0.1mol/L CeSO4</b>	<b>8/10/2012 10:44:00 AM</b>
<b>Start time:</b>	<b>8/10/2012 11:11:10 AM</b>		

---

**Titer**

Titrant	Ce(SO4)2 c = 0.1 mol/L
Titer	0.99808

- 
- (1) Modified
  - (2) Excluded
  - (3) Outside limits
  - (4) Resource expired
  - (5) srel above max srel
  - (6) srel above max srel for multiple determination
  - (7) Value outside limits, not saved in setup
  - (8) Sample data outside limits
  - (9) Standard evaluation used
  - (10) Result from buffer

**Created:**                    **Development Administrator (admin), 8/10/2012 11:17:08 AM**